# **ELBER**

## M E D I C A L

MANUAL DE INSTRUÇÕES E GARANTIA  
INSTRUCTION AND WARRANTY MANUAL  
MANUAL DE INSTRUCCIONES Y GARANTÍA

# CAMARA CONSERVADORA ELBER PARA TRANSPORTE

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14	GARANTIA .....	16
15	CONSERVATOR CHARACTERISTICS .....	17
15.1	Packaging.....	18
15.2	Transportation .....	18
15.3	Storage .....	18
15.4	Handling.....	18
16	CONSERVATIVE INSTALLATION AND FIXATION .....	18
17	ELECTRICAL INSTALLATION .....	19
18	OPERATION .....	21
19	SOUND ALARM.....	21
20	BATTERY.....	22
21	BEST CONSERVATOR USE AND RECOMMENDATIONS.....	22
22	DEGREE AND CLEANING .....	23
23	MAINTENANCE AND PREVENTIVE VERIFICATION.....	23
24	EMERGENCY SYSTEM (see optional availability) .....	24
25	AUTO ELECTRONIC DIAGNOSTICS (Only applicable for conservatives with 12/24V compressor).....	25
26	PROBLEMS AND PROBABLE CAUSES.....	25
27	DISPOSAL OF THE PRODUCT .....	26
28	WARRANTY .....	26
29	CARACTERÍSTICAS DE LA CONSERVADORA .....	27

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## **14 GARANTIA**

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## **Dear USER!**

You have just purchased a preservative for transporting immunobiologicals, blood products, medicines, thermolabile and the like. This product has been developed with ELBER technology and functionality. We are grateful for your choice.

It is essential to read carefully the instructions contained in this manual.

Following these guidelines, you can better understand the operation of the product, taking advantage of its benefits and extending its durability. The instruction manual must always be accessible near the product so that any doubt can be remedied immediately.

## **15 CONSERVATOR CHARACTERISTICS**

ELBER transport preservatives have been developed to store immunobiological, blood derivatives, medicaments, thermolabile and the like at a stable, homogeneous and uninterrupted temperature. It is indicated for automotive installation (cars, vans, buses, etc).

Operates directly on the vehicle battery (12 or 24VDC).

The conservative ELBERs are designed to ensure that the internal chamber temperature is always stable and homogeneous. Equipped with safety system, the conservative Elber remain working even with the vehicle off (see option availability). The conservative Elber has an electronic microprocessor controller that controls the internal temperature of the chamber (cooling as needed) and informs the user of all events through audible and visual alarms and electronic records.

For better internal comfort of the vehicle, this conservator is mounted with UCR (remote compressor unit).

**REMOTE COMPRESSOR UNIT:** The compressor (engine) of this conservator is installed separately, inside the car. The refrigeration cabinet is connected to the remote compressor unit, using an Elber flexible multi-tube hose. This hose is coupled to the compressor unit, through quick couplings.

## 15.1 Packaging

The packaging of the product you received was designed to provide security in the conservative transport. When unpacking, check the integrity of the product. The conservative ELBER are packed with bubble wrap and cardboard. Some products have packaging in MDF or wooden pallet. Each product is individually packaged.

After unpacking the product, all packing material should be sent for recycling or reuse.

## 15.2 Transportation

Conservative ELBERs must be transported in an upright position, maintaining the packaging conditions according to the label.

## 15.3 Storage

ELBER preservatives should be stored in clean, dry, sun-protected and excessive humidity environments.

**The preservative should not be stored for more than three months without being used.** If it is not used, the battery will discharge (see option availability).

## 15.4 Handling

ELBER preservatives should only be handled by trained professionals who fully meet the specifications in the product instruction manual.

## 16 CONSERVATIVE INSTALLATION AND FIXATION

The retainer should be fixed to the vehicle by means of bolts, fixed in resistant parts of the car. If necessary, the external cabinet can be used for fixing. Give preference to the corners, which have greater resistance.

**ATTENTION:** Do not insert drills or screws longer than 15 mm as this may affect the refrigerant piping.

For products with a specific fastening system, see instruction on the technical specifications sheet.


For conservatories fitted with a water drain, keep the drain cap closed and only open during cleaning and draining of the product.

The conservative Elber should be installed preferably protected from the sun and in a place where there is air circulation. Very small rooms should have openings for air exchange.

We do not recommend installing the conservator near heat sources (heaters and greenhouses). In the ventilation of the compressor unit (UCR) there must be space of at least 20 centimeters from the wall or any obstacle, to allow ventilation of the compressor unit and the electronic system.

It can be built into furniture, as long as its remote compressor unit is not obstructed, preserving the opening of the free ventilation grids for air circulation. This is necessary so that there is heat exchange with the environment and free space to remove the lid that gives access to the refrigerator compressor.

**Ventilation of the compressor unit can not be obstructed. The heat exchange of the compressor unit is indispensable for the proper functioning of the product.**

The ventilator flaps of the preservative are indicated on the product data sheet.

Observe the identification of air inlet and outlet. The air enters through an opening and leaves through another. This flow is promoted by the fan that is installed inside the compartment.

## **17 ELECTRICAL INSTALLATION**

The conservative Elber are designed to be powered by the vehicle's battery. Its voltage is factory set at 12, 24 or 12/24 VDC.

**Check the supply voltage on the product identification label.**


**The polarity indicated on the harness must be respected.** Do not consider the coloring of the wires to guide the polarity. In the refrigerator the blue color is positive and in the vehicle the blue color may be negative. Consider the (+) and (-) identification washers on the cables.

Reversing the polarity may cause the electronic unit or the protective fuse to burn.

To connect the conservative Elber, simply plug the conservator plug into the power plug. For refrigerators that do not have a plug, simply patch the cables.

For products with specific electrical connection, check the electrical connection diagram on the technical specifications sheet of the product.

For other products, it is recommended that you follow the wiring diagram below, which helps to maintain the conservator / auxiliary battery charge and does not discharge the car battery. Some conservatives have a battery charger, which can be operated with the vehicle stopped (see supply voltage on the identification label).

To install the conservator, use the electric conductor of the battery to power the refrigerator. Dimension the electric conductor (wire) according to the distance between the conservator and the battery according to the following table:

	<i>BITOLA</i>		12 V	24 V
CONDUCTOR	2,5 mm <sup>2</sup>	For distance to	2,5 m	5 m
CONDUCTOR	4,0 mm <sup>2</sup>	For distance to	4 m	8 m
CONDUCTOR	6,0 mm <sup>2</sup>	For distance to	6 m	12 m
CONDUCTOR	10,0 mm <sup>2</sup>	For distance to	10 m	20 m

The above table is indicated for products with 1 compressor, for products with more than 1 compressor, consult AT Elber to obtain the correct specification.

## 18 OPERATION

This preservative is prepared for the preservation of immunobiological, blood derivatives, medicaments, thermolabile and the like.

It guarantees a constant temperature at any time of the year and works even when there is a lack of electricity (see optional availability). The temperature is adjusted and controlled automatically by the system and there is no need for intervention even on cold or hot days.

**Ensuring uniform temperature** is made by the digital controller that maintains what has been programmed. The controller is programmed to keep the chamber at **+ 2.0 ° C to + 8.0 ° C** on both cold and hot days. At any time of the year, winter or summer, the conservative will adjust automatically.

The temperature will always remain in a defined range, ensuring the best properties of the preserved material.

## 19 SOUND ALARM

For convenience, the conservator has an audible alarm that is triggered when the temperature is outside the programmed limits or other anomalies. This can occur if the door is too long open or due to any malfunction. So in a timely manner some action may be taken so that the products in conservation are not lost. The alarm will only be activated after the preset time with a temperature outside the programmed time.

**Delay alarm function:** When the conservator is turned on, the alarm will sound. To suspend the alarm for the preset time, press the muting key (ELV5 and touch screen) or any key (ELV4) of the digital controller. Repeat as many times as necessary until the conservator reaches the programmed range, where it will no longer fire. To change the time between shots, see regular time function to trigger the alarm.

## 20 BATTERY

If the product has an emergency system, it is important to follow the recommendations below:

- For conservatives with emergency system and are purchased without the battery, always use new automotive battery.
- When using unsealed battery, check distilled water level quarterly.
- Verificar o funcionamento da bateria. Caso não esteja funcionando corretamente, ocasionará o mau funcionamento da conservadora.

## 21 BEST CONSERVATOR USE AND RECOMMENDATIONS

- Avoid leaving the door open;
- In use, always keep products inside the conservator, even if they are bottled with water. This load of products helps keep the temperature inside the conservator stable;
- Before putting the vaccine in the conservator check that the operation is correct and the air circulation fan in the inner bottom is working;
- Avoid leaving the conservator door open for long or open it constantly;
- Keep the conservator running for at least 48 hours before placing products in the conservator to verify that the general operation is correct, if the internal fan is circulating the air and if the temperature is within the set parameters;
- In use, always keep products inside the conservator, even if they are bottled with water. This load of products inside helps to maintain stable the temperature inside the conservative;
- The vents of internal ventilation can not be obstructed, constantly check if stored products are obstructing the fins and, therefore, the internal ventilation of the product;
- The shelves or drawers can not be lined with plastics, towels or the like, which compromises the internal ventilation of the product;
- On days with high humidity, there may be condensation of air moisture (sweat) on the door or on the outer edges of the conservator (near gasket). In places with room temperature controlled by air conditioning, this characteristic is almost extinct;
- Para conservadoras equipadas com sistema de emergência, não realizar battery runtime tests without Elber technical support, on-site or on-line. Only technicians can guarantee the recording of results, disable the messaging system and ensure that the battery is not damaged (*see optional availability*);
- Only trained personnel should operate the product.

## 22 DEGREE AND CLEANING

To clean the conservator, first turn off at the GENERAL SWITCH. Do not simply unplug it as it will continue to run through the battery (see option availability).

Use damp cloths with water and mild detergent. Do not use splashes or water hoses. Excess water may compromise product performance.

Do not use chemicals (alcohol, solvents, etc.) or abrasive products and hard sponges. Non-recommended procedures may stain and damage the surfaces and coating of the refrigerator.

During cleaning, remove all materials stored in the conservator and place them in another equipment or box with cryogenic material.

**Important!** After a period of inactivity of the conservator (for cleaning, for example), when the power is turned on, the conservator can reach temperatures below that programmed in the electronic controller. Therefore, whenever you turn off the conservator, after reconnecting, wait at least 2 hours to put load on the conservator again. After 2 hours, clear the controller data and start normal use.

## 23 MAINTENANCE AND PREVENTIVE VERIFICATION

- Monthly check that the gasket of the door gasket is correct (there must be no gaps between the gasket and the cabinet);
- Monthly check the level of glycerol present in the bottle, located inside the conservator (variable position depending on the model, some models do not have the sensor immersed in glycerol). If the sensor is out of solution, complete the level;
- Monthly visually verify the presence of cracks or broken parts in the conservator's structure;
- Monthly check the tightening of the door hinge bolts manually;
- Monthly visually check the integrity of door locks. Perform practical test of the operation of the latches;
- Monthly visually verify the integrity of the cushions and the conservative base of support (when it exists);
- Monthly visually check if the electrical installation is correct (presence of oxidation in plugs, relays and components in general, tightening of the terminals and plugs, state of electrical cables, etc);

- For conservators equipped with an emergency system, check battery charge every six months (see optional availability);
- Conservative condenser should be cleaned annually. The location of the condenser varies for each model, but is always close to the compressor;
- For conservatives equipped with cryogenic blocks, cryogenic blocks need to be replaced after expiry date. The expiry date is stated on the carton of the cryogenic block (see optional availability);
- Change the evaporator and condenser fans of the conservator every 4 years;
- For conservatives equipped with an emergency system, change the conservator's battery every 4 years (see optional availability).

**NOTE: Monitoring systems that depend on the transmission and/or storage of data by wifi, gps, sms, WhatsApp or satellite are subject to changes in technologies, licenses and possible related costs, to monitor maintenance.**

## **24 EMERGENCY SYSTEM (see optional availability)**

When equipped with an emergency system, in the absence of vehicle power (12 or 24VDC), the Elber preservatives continue to operate (see the technical specifications sheet of the product), keeping all their electrical and electronic functions in operation.

The activation of the emergency system is automatic. The autonomy of the emergency system varies for each model, can be increased or decreased, according to the need of the client.

When equipped, the Elber Conservatives have sealed, maintenance-free batteries with an automatic charger.

Battery recharging time varies between 10 and 20 hours, and this time may be longer if the battery bank is too large.

**Important:** The autonomy of the conservator depends on many factors, such as: battery capacity used, ambient temperature, conservative load, number of times the door is opened, among others. Use as reference a 110 Ah automotive battery, which should have an autonomy of approximately 10 hours (considering ambient temperature  $22 \pm 2^{\circ}\text{C}$ , few door openings, new and conservative automotive battery with maximum load already at standard temperature). The higher the battery amperage, the greater the autonomy of the conservator.

The emergency system is optional, check the availability of the optional. Some models can be sold with the emergency system, however, without the battery.

## 25 AUTO ELECTRONIC DIAGNOSTICS (Only applicable for conservatives with 12/24V compressor)

LED	Blinking	Diagnosis
Blinking	1 X	Battery cut-off (voltage not specified)
Blinking	2 X	Many attempts to start the compressor or short fan (current above 0.5A)
Blinking	3 X	Compressor did not break or over pressure. Wait 4 minutes for it to automatically reconnect.
Blinking	4 X	Minimum compressor rotation error
Blinking	5 X	Ambient temperature too high
Blinking	6 X	Thermostat is shorted or the electronic unit is incorrectly connected.

## 26 PROBLEMS AND PROBABLE CAUSES

FAILURE	SOLUTION
Conservative does not turn on	<ul style="list-style-type: none"> <li>- Check that the supply voltage is the same as that of the conservator;</li> <li>- Check for power in the power supply;</li> <li>- Check that the main switch is switched on;</li> <li>- Check that the power plug is securely connected;</li> <li>- Check that the fuse is blown;</li> <li>- Check that the set temperature in the controller is correct;</li> <li>- Check if there is any warning on the digital controller.</li> </ul>
Conservative does not cool	<ul style="list-style-type: none"> <li>- Check that the compressor is switched on (place your hand on the compressor and check for vibration and slight heating);</li> <li>- Check that the internal ventilation is working;</li> <li>- Check the door closing (door gasket must not have cracks in relation to the cabinet);</li> <li>- Check if there is any warning on the digital controller.</li> </ul>
Conservative alarm by high temperature	<ul style="list-style-type: none"> <li>- Make sure that a large load of product has been charged within 2 hours (the conservator can take up to the optimum temperature when a large quantity of product is placed at one time);</li> <li>- Check door closing (door gasket must not have cracks in relation to cabinet).</li> <li>- Check if there is any warning on the digital controller.</li> </ul>
Conservative internal temperature varying	<ul style="list-style-type: none"> <li>- Check that the internal fans are working;</li> <li>- Check that the internal ventilation flaps are free (they can not be obstructed);</li> </ul>

## **27 DISPOSAL OF THE PRODUCT**

Disused or abandoned refrigerators can be dangerous, especially for children, as they can get stuck inside, running the risk of shortness of breath. Before discarding your old refrigerator, cut off the power cord and remove the door (s). Leave the shelves in place so the children can not enter.

If the conservator's battery is discarded, route the discarded battery to a nearby pickup point.

In case of product disposal, all materials should be sent for recycling.

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