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REFRIGERADOR CFMB

ELBER

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34	PREVENTIVE MAINTENANCE AND VERIFICATION	38
35	PROBLEMS AND PROBABLE CAUSES.....	39
36	PRODUCT DISPOSAL.....	39
37	WARRANTY	39
38	PROCEDURE FOR HANDLING THE 9-NUMBER ALARM DIALER (<i>SEE OPTIONAL AVAILABILITY</i>).....	40
38.1	Programming a Phone Number	40
38.2	Programming a phone number with pause 2 secs. (for PBX calls)	40
38.3	Deleting a phone number.....	40
38.4	Testing the phone numbers	41
38.5	Dialing process.....	41
39	SURGE PROTECTION (DPS).....	41
39.1	How DPS works.....	42
39.2	Connection and attachment of the ELBER DPS kit	42
40	TERMO DE GARANTIA.....	44

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DEAR USER!

You have just acquired a dual-compartment equipment designed for the storage of immunobiologicals, blood products, medicines, thermolabiles and the like. This product was developed with ELBER technology and functionality. We are grateful for your choice.

It is essential to carefully read the instructions contained in this manual.

By following these guidelines, you will be able to better understand how the product works, enjoying its benefits and extending its durability. The instruction manual should always be accessible near the product, so that any doubts can be answered immediately.

ANVISA Registration: 80698750001, 80698759001, 80698750002, 80698750003, 80698750004

ANVISA Authorization (AFE): 8069875.

27 FREEZER FEATURES

ELBER products have been developed to store at a stable and homogeneous temperature.

It works through electrical power 110-127 VAC or 220-240 VAC.

ELBER products have been designed to ensure that the internal temperature of the chamber is always stable and homogeneous. These equipment have a microprocessor electronic controller that controls the internal temperature (cooling as needed) and informs the user of all events through audible and visual alarms and electronic records.

ELBER products can be equipped with optional amusements, see the technical specification sheet, which accompanies your product, all the features of the product.

27.1 Packing

The packaging of the product you received is designed to provide safety in the transport of the equipment. When unpacking, check the integrity of the product. ELBER products are packed with bubble wrap and cardboard. Some products have MDF or wooden bed frame packaging. Each product is individually packaged.

After unpacking the product, any and all packaging material must be sent for recycling or subsequent reuse.

27.2 Transport

ELBER equipment must be transported in an upright position, maintaining packaging conditions according to the label.

27.3 Storage

ELBER equipment should be stored in clean, dry environments, protected from the sun and excessive humidity.

27.4 Handling

ELBER equipment should only be handled by trained professionals who fully meet the specifications in the product's instruction manual.

28 INSTALLATION AND FIXING


The equipment must be positioned on a flat and level floor. The feet have casters that allow it to be positioned easily. To move the equipment, the casters must be unlocked and locked when the equipment is positioned. Moving it with the casters locked may damage the casters and the floor.

It should preferably be installed protected from the sun and in a place where there is air circulation. Very small rooms should have openings for air exchange.

We do not recommend installing the equipment near heat sources (heaters and stoves).

In the ventilation grilles of the equipment, there must be a space of at least 20 centimeters from the wall or any obstacle, to allow ventilation of the compressor unit and the electronic system.

It is recommended that the equipment be installed in a room with air conditioning. As the equipment works at extremely low temperatures, air conditioning helps the equipment to maintain its good performance and reduces the relative humidity of the air, which prevents external condensation events in the equipment.

29 WIRING

The equipment is prepared to work connected to 110-127 VAC or 220-240 VAC power.

Before switching on, check that the voltage indicated on the product coincides with the mains voltage.

The power cord must be plugged into a permanent wall or countertop socket outlet in accordance with current regulations. Do not use extension cords or adapters.

Observe the steps below to turn on the equipment correctly.

1st - The voltage 110-127 VAC or 220-240 VAC is set at the factory,

2nd - Plug the plug into the socket - the equipment will not work yet, as the main switch is still to be turned on;

3rd - Turn on the equipment's main switch;

4th - When turning on the equipment, the **Audible Alarm** will sound until the temperature inside the chamber is within the ranges established in the digital controller.

You can cancel the audible alarm for a predetermined time by pressing the mute key. When the alarm is sounding and is deactivated, it will sound again after the programmed time if the temperature remains outside the programmed range. The alarm will sound until the conservator reaches the normal operating range (on average 2 hours without load).

To turn off the conservator, when equipped with an emergency system, it is necessary to turn off the main switch. It is not enough to remove the plug from the socket, as the conservator will continue to work on the battery (*see availability of the optional*).

The use of the SPW 275-45 surge protection device (SPD) is indispensable. This device protects the conservator against surges and over voltages, preventing it from burning.

Attention!

This conservator must be connected to a ground to avoid damage to the product and the user. Ground wire connection is required at any voltage (110-124 VAC or 220-240 VAC). The lack of grounding and use of the DPS device compromises the coverage of the conservator's warranty. To perform the grounding and electrical installation, as well as the SPD, follow the standards of ABNT-NBR 5410, or call a responsible professional.

30 OPERATION

Once powered on, the Elber equipment takes about 2 hours to reach the programmed temperature range (*). The controller is already factory adjusted and does not require any type of adjustment, just turn on the product and wait to reach the working temperature range.

The **guarantee of a uniform temperature** is made by the digital controller, which maintains what has been programmed. The controller is programmed to keep the camera at a temperature between **+2.0°C to +8.0°C** for the conservator compartment and **-30.0°C to -20.0°C** for the freezer compartment. The temperature will always remain in a defined range, ensuring the best material properties inside the equipment.

() No load, ambient temperature of 22°C and no door openings.*

31 SOUND ALARM

For added convenience, the equipment has an audible alarm that is triggered when the temperature is outside the programmed limits or other anomalies. This can occur if the door is left open for a long time or due to some malfunction. Thus, in a timely manner, some measure can be taken so that the products in conservation are not lost. The alarm will activate only after the predetermined time with temperature outside the programmed range.

Alarm delay function: When the equipment is turned on, the alarm will sound. To suspend the alarm for the predetermined time, press the digital controller's mute key. Repeat the operation as many times as necessary until the equipment reaches the programmed range, where it will no longer fire. To change the time between shots, see Regular time function to trigger the alarm.

32 BEST USE OF EQUIPMENT AND RECOMMENDATIONS

- Keep the conservator in an air-conditioned environment (+22 °C).
- In use, always keep products inside the conservator, even if they are bottles with water. This load of products helps to keep the temperature inside the conservator stable.
- Before placing the vaccine in the conservator, check that the operation is correct and the air circulation fan at the bottom inside is working.
- Avoid leaving the conservator door open for a long time or opening it constantly;

- Keep the conservator in operation for at least 48 hours before placing products in the conservator, to check that the general operation is correct, that the internal fan is circulating air and that the temperature is within the regulated parameters;
- The internal ventilation fins cannot be obstructed, constantly check if stored products are obstructing the fins and, consequently, the internal ventilation of the product;
- Shelves or drawers cannot be lined with plastics, towels or the like, which compromises the internal ventilation of the product;
- On days with a lot of humidity, there may be condensation of air humidity (sweat) on the door or on the outer edges of the conservator (near the gasket). In places with ambient temperature controlled by air conditioning, this characteristic is practically extinct.
- Only trained persons should operate the product.

33 DEFROSTING AND CLEANING

The conservative compartment has an automatic defrosting system, there is no need for user intervention, that is, the product adjusts to the use and makes automatic defrosts as needed.

The freezer compartment has manual defrosting, that is, it is necessary to monitor the thickness of ice on the internal walls of the product. When this wall is equal to or greater than 10 mm, it is recommended to defrost the product. To defrost, it is necessary to remove all the stored material, turn off the product (turn off the main switch and unplug it) and leave the door open. The ice accumulated on the walls will slowly melt. It is also necessary to remove excess water at the bottom of the product as the thaw happens, using dry cloths or absorbent sponges.

Attention: It is not enough to unplug the product from the socket, as the controller will continue to operate on battery.

For cleaning use damp cloths with water and neutral detergent. Do not use water nozzles or hoses, excess water may compromise the operation of the product.

Do not use chemical products (alcohol, solvents, etc.) or abrasive products and hard sponges. Procedures not recommended may stain and damage the surfaces and coating of the equipment.

During cleaning, remove all stored materials and place them in another equipment or cooler with cryogenic material.

34 PREVENTIVE MAINTENANCE AND VERIFICATION

- Monthly check if the sealing of the door gaskets is correct (there can be no gaps between the gasket and the cabinet);
- Monthly check the level of glycerol present in the bottle, located inside the equipment (variable position depending on the model, some models do not have the sensor immersed in glycerol). If the sensor is out of solution, complete the level;
- Annually the condenser of the equipment must be cleaned. The location of the condenser varies for each model, but it is always close to the compressor;
- From 12 months on, it is important to check for soot or dirt that may be clogging the condenser;
- After 03 years of use of the equipment, it is recommended to test the battery life of the emergency system (test consists of removing the equipment from the outlet for 2 hours, it must be monitored, to see if there is an oscillation in the internal temperature of the equipment) Periodicity – every 6 months;
- After 05 years of use of the equipment, it is recommended to replace the internal fans and condenser fan of the chamber;

According to the special characteristics of the equipment and/or customer requirements, systematic preventive maintenance can be adapted;

Monitoring systems that depend on transmission and/or storage of data by Wifi, GPS, SMS, WhatsApp, or via satellite are subject to changes in technologies, licensing and possible related costs, to maintain monitoring.

35 PROBLEMS AND PROBABLE CAUSES

FAULT	SOLUTION
The equipment does not turn on	- Check if the voltage of the commercial electrical network is the same as the equipment; - Check if there is commercial electricity; - Check if the main switch is turned on; - Check if the plug is plugged into the socket; - Check if there is any warning on the digital controller.
The equipment does not freeze	- Check if the compressor is on (place your hand on the compressor and check for vibration and slight heating); - Check if the internal ventilation is working (when applicable); - Check the closing of the door (door gasket cannot have gaps in relation to the cabinet); - Check if there is any warning on the digital controller.
Equipment is alarming by high temperature	- Check if a large load of product has been charged for less than 2 hours (the equipment can take time to reach the ideal temperature when a large amount of product is placed at once); - Check if the door is being opened too often and/or for too long periods; - Check the closing of the door (door gasket must not have gaps in relation to the cabinet). - Check if there is any warning on the digital controller.
Internal temperature of the equipment varying	- Check if the internal fan is working; - Check that the internal ventilation fins are free (they cannot be obstructed); - Check if the door is being opened too often and/or for too long periods; - Check if there is any warning on the digital controller.

36 PRODUCT DISPOSAL

Unused or abandoned refrigerators can be dangerous, especially for children, as they can get trapped inside, running the risk of shortness of breath. Before disposing of your old refrigerator, cut the power cord and remove the door(s). Leave the shelves in place so that children cannot enter.

In case of disposal of the battery of the equipment, send the discarded battery to the nearest collection point.

In case of product disposal, all materials must be sent for recycling.

37 WARRANTY

Equipment warranty is 1 year.

38 PROCEDURE FOR HANDLING THE 9-NUMBER ALARM DIALER *(SEE OPTIONAL AVAILABILITY).*

The tele alarm works through a telephone line. It makes calls to the registered numbers when the temperature is out of the programmed range, if the battery voltage is low, over mains operation and open door. For the dialer to work, you need to connect a telephone cable.

For more information check the specification of your equipment.

38.1 Programming a Phone Number

Press # + Installer Password + Index (01 to 09) + Phone Number (up to 22 digits) + # + #.

38.2 Programming a phone number with pause 2 secs. (for PBX calls)

Press # + Installer Password + Table of Contents (01 to 09) + 0 2 second pause + E + Phone Number + # + #.

38.3 Deleting a phone number

Press # + Installer Password + Number Position you want to delete (Index 01 to 09) + * + # + #.

Another way to delete a programmed number is to access the desired index and omit the value of the phone number by pressing #. Be careful when accessing the phone indexes (01 to 09) accidentally, as you will need to re-enter their value.

38.4 Testing the phone numbers

With the telephone line properly connected and the dialer in user mode (out of programming), press F + Desired Memory Position (1 to 9 – without the digit 0).

The dialer will place a call to the stored number, playing a siren sound, on the phone, for 60 seconds. During this time, the LINE LED will remain on. To stop a test in progress, press any key.

38.5 Dialing process

The dialer goes through a call cycle that begins with the phone registered in the first memory position (index 01) and ends with the phone registered in the last memory position (index 09), performing up to 2 dialing cycles.

Ex: Start of the first cycle / MEM1... MEM2... MEM3... MEM4... MEM5... MEM6... MEM7... MEM8... MEM9... / End of the first cycle / MEM1... MEM2... MEM3... MEM4... MEM5... MEM6... MEM7... MEM8... MEM9... / End of the second cycle and dialing process.

THE OTHER FUNCTIONS OF THE TELE ALARM ARE ALREADY PROGRAMMED AT THE FACTORY, AND IT IS NOT NECESSARY TO CHANGE.

39 SURGE PROTECTION (DPS).

The Elber SPD is a surge protection device (SPD). It is an equipment that aims to detect overvoltages in the electrical network and divert surge currents, it is an indispensable device for the use of Elber products. It can be installed in the circuit breaker box present on site or through the ELBER DPS Kit.

Maximum device discharge current: 45 kA

39.1 How DPS works


The Elber SPD is connected to the power grid and responsible for surge protection, when subjected to an "Electrical Surge" the SPD will stabilize the electrical current that the equipment will receive.

ATTENTION: The Surge Protector when subjected to a very large power variation can be damaged, it will signal through an LED that will light up, and the replacement of the protector will be necessary.

ATTENTION: to perform the installation of the SPD follow the standards of ABNT-NBR 5410:2004, section DPS, or call a responsible and qualified professional.

39.2 Connection and attachment of the ELBER DPS kit

To connect the DPS Kit, it is essential to turn off the electric current!


ELBER DPS KIT LINK:

The Elber DPS kit consists of a Surge Protector, a female socket for connecting the equipment to be connected and a box.

The Elber DPS kit must be connected directly to the mains, respecting all required safety standards. It is essential to connect the "Terra" to the kit, as it can provide all the security proposed to your equipment.


The kit comes with a cable for connecting the Phase, Neutral and Earth, which will come out from the bottom of the DPS kit.

ATTENTION: It is essential to correctly insulate the wires, either by insulating tape or heat shrink tubing.

DPS KIT ATTACHMENT:

After connecting the SPD to the mains, and still with the power off, remove the top cover of the Kit box, if necessary, also remove the Surge Protector, for better access to the fixing points located at the bottom of the box.

To attach the DPS Kit box, simply position it to the desired location and insert the screws through the attachment points.

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This product is insured against manufacturing and/or functional defects for 12 months after delivery to the customer.

All parts and components are covered by the warranty, except when damaged by accidents, recklessness or failure to observe the installation instructions for use of the equipment.

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